





ARTHROSCOPY IMPLANTS & INSTRUMENTS

Cannulated Interference Screw (ACL Screw)

Cannulated Interference Screw (ACL Screw) (Sterile)

(Sterile)

Code No.	Diameter	Code No.	Diameter
T315A	7.0mm	ST315A	7.0mm
T315B	8.0mm	ST315B	8.0mm
T315C	9.0mm	ST315C	9.0mm
T315D	10.0mm	ST315D	10.0mm



LENGTH

20mm to 35mm (5mm Increment)

MATERIAL

Titanium

T316A Herbert Screw 2.4mm Cannulated

ST316A Herbert Screw 2.4mm Cannulated (Sterile)



LENGTH

14mm to 30mm (2mm Increment)

MATERIAL

Titanium

T316 Herbert Screw 3.5mm Cannulated

ST316 Herbert Screw 3.5mm Cannulated (Sterile)



LENGTH

14mm to 40mm (2mm Increment)

MATERIAL

Titanium

T2001 Headless Screw 2.5mm Cannulated



LENGTH

08mm to 30mm (2mm Increment)

MATERIAL

Titanium

T2002 Headless Screw 3.0mm Cannulated




LENGTH

16mm to 30mm (2mm Increment)
35mm to 45mm (5mm Increment)

MATERIAL

Titanium

T2003 Headless Screw 3.5mm Cannulated




LENGTH

16mm to 30mm (2mm Increment)
35mm to 50mm (5mm Increment)

MATERIAL

Titanium

T2004 Headless Screw 4.0mm Cannulated




LENGTH

16mm to 30mm (2mm Increment)
35mm to 50mm (5mm Increment)

MATERIAL

Titanium

T2005 Headless Screw 5.0mm Cannulated



LENGTH

30mm to 80mm (5mm Increment)

MATERIAL

Titanium

T317.01 Endo Buttons

ST317.01 Endo Buttons (Sterile)



LENGTH

15mm to 30mm (5mm Increment)

MATERIAL

Titanium

ST317 Endo Buttons with Loop (Sterile)




LENGTH

15mm to 30mm (5mm Increment)

MATERIAL

Titanium

ST317.02 Endo Buttons Adjustable (Sterile)



MATERIAL

Titanium

T318.01 Suture Disc

ST318.01 Suture Disc (Sterile)



MATERIAL

Titanium

T323 Anchor Screws

ST323 Anchor Screws (Sterile)



DIAMETER

3.0mm, 5.0mm

MATERIAL

Titanium

ST323 Suture Anchor Screws (Sterile)



DIAMETER

3.0mm, 5.0mm

MATERIAL

Titanium

ST323A Suture Anchor Screws with Needle (Sterile)



DIAMETER

3.0mm, 5.0mm

MATERIAL

Titanium

T340 Eight Plate



LENGTH

12mm to 18mm (2mm Increment)

MATERIAL

Titanium

Implants

T340A Eight Plate Z Angle



LENGTH

12mm to 18mm (2mm Increment)

MATERIAL

Titanium

T384 Spiked Washer



SIZE

4.0mm (5 Spikes), 4.5mm (7 Spikes), 6.5mm (9 Spikes)

MATERIAL

Titanium

T385 Spiked Ligament Staple



LENGTH

6.0mm, 8.0mm, 10.0mm

MATERIAL

Titanium

673.001 Main zig



673.002 Acl elbow aimer



673.004 Pcl Femoral



673.005 Pcl Tibia



673.006 Angle Bullet




673.007 Graft passing 2.4mm drill bit




661.005 Beth Pin Drill tip



663 ACL Tibia Straight Cannulated Reamer
6mm, 7mm, 8mm, 9mm, 10mm, 11mm



- 662** ACL Femoral Reamer Cannulated Flower Tip
6mm, 7mm, 8mm, 9mm, 10mm, 11mm



- 661.002** Endoscopic Cannulated Reamer



- 673.008** Depth Gauge for ACL/PCL reconstruction



- 666.050** Femoral Aimer With Offset 5mm

- 666.070** Femoral Aimer With Offset 7mm



- 661.078** Tunnel Plug, 7-8mm

- 661.089** Tunnel Plug, 8-9mm

- 661.910** Tunnel Plug, 9-10mm



- 661.070** Tunnel Dilator 7mm

- 661.080** Tunnel Dilator 8mm

- 661.090** Tunnel Dilator 9mm



- 659.035** Screw Driver For Interference Screw



- 673.009** Impactor for spiked ligament staple



- 673.010** Arthroscopic Probe




- 674.030** Microfracture Awl 30'

- 674.045** Microfracture Awl 45'

- 674.090** Microfracture Awl 90'



661.003 Graft Master Board, Suturewise with Tensiometer



668.001 Graft Sizing Block



665.060 Tendon Stripper Close 6mm


665.070 Tendon Stripper Close 7mm



674.001 Straight Rasp



674.002 Convex Rasp



674.005 Femoral PCL Elevator



674.006 T-Handle for Tunnel Dilators



SAGC.45 Knee Arthroscopy Instrument Set



Code No.	Description	Qty.
673.001	Main zig	1
673.002	Acl elbow aimer	1
673.004	Pcl Femoral	1
673.005	Pcl Tibia	1
673.006	Angle Bullet	1
673.007	Graft passing 2.4mmdrillbit	1
661.005	Beth Pin Drill tip	1
663.060	ACL Tibia Straight Cannulated Reamer 6mm	1
663.070	ACL Tibia Straight Cannulated Reamer 7mm	1
663.080	ACL Tibia Straight Cannulated Reamer 8mm	1
663.090	ACL Tibia Straight Cannulated Reamer 9mm	1
663.100	ACL Tibia Straight Cannulated Reamer 10mm	1
663.110	ACL Tibia Straight Cannulated Reamer 11mm	1
662.060	ACL Femoral Reamer Cannulated Flower Tip 6mm	1
662.070	ACL Femoral Reamer Cannulated Flower Tip 7mm	1
662.080	ACL Femoral Reamer Cannulated Flower Tip 8mm	1
662.090	ACL Femoral Reamer Cannulated Flower Tip 9mm	1
662.100	ACL Femoral Reamer Cannulated Flower Tip 10mm	1
662.110	ACL Femoral Reamer Cannulated Flower Tip 11mm	1
661.002	Endoscopic Cannulated Reamer	1
673.008	Depth Gauge for ACL/PCL reconstruction	1
666.050	Femoral Aimer With Offset 5mm	1

Code No.	Description	Qty.
666.070	Femoral Aimer With Offset 7mm	1
661.078	Tunnel Plug, 7-8 mm	1
661.089	Tunnel Plug, 8-9 mm	1
661.910	Tunnel Plug, 9-10 mm	1
661.070	Tunnel Dilator, 7mm	1
661.080	Tunnel Dilator, 8mm	1
661.090	Tunnel Dilator, 9mm	1
659.035	Screw Driver For Interference Screw	1
673.009	Impactor for spiked ligament staple	1
673.010	Arthroscopic Probe	1
674.030	Microfracture Awl, 30°	1
674.045	Microfracture Awl, 45°	1
674.090	Microfracture Awl, 90°	1
661.003	Graft Master Board, Suturewise with Tensiometer,	1
668.001	Graft Sizing Block	1
665.060	Tendon Stripper Close 6mm	1
664.070	Tendon Stripper Open 7mm	1
674.001	Straight Rasp	1
674.002	Convex Rasp	1
674.005	Femoral PCL Elevator	1
674.006	T-Handle for Tunnel Dilators	1

669.001 Jig for Herbert Screw



658.024 Herbert Screw Tap Cannulated 2.4mm

658.035 Herbert Screw Tap Cannulated 3.5mm

658.045 Herbert Screw Tap Cannulated 4.5mm



659.025 Screw Driver Cannulated 2.5mm Tip



659.020 Screw Driver Cannulated 2.0mm Tip



657.020 Drill Bit Cannulated 2.0mm

657.027 Drill Bit Cannulated 2.7mm

657.032 Drill Bit Cannulated 3.2mm



671.100 Guide Wire 1.0mm x 225mm

671.125 Guide Wire 1.25mm x 225mm



670.125 Direct Measuring Device for 1.25mm Guide Wire



SAGC.52 2.4/3.5mm Herbert Screw and Instrument Set



Instruments

Code No.	Description	Qty.
669.001	Jig for Herbert Screw	1
658.035	Herbert Screw Tap Cannulated 3.5mm for upper thread of Herbert Screw	1
658.045	Herbert Screw Tap Cannulated 4.5mm for Lower thread of Herbert Screw	1
658.024	Herbert Screw Tap Cannulated 2.4mm for Lower thread of Herbert Screw	1
659.025	Screw Driver Cannulated 2.5mm Tip	1
659.020	Screw Driver Cannulated 2.0mm Tip	1
657.027	Drill Bit Cannulated 2.7mm	1
657.032	Drill Bit Cannulated 3.2mm	1
657.020	Drill Bit Cannulated 2.0mm	1
670.125	Direct Measuring Device for 1.25mm Guide Wire	1
671.125	Guide Wire 1.25mm x 225mm	3
671.100	Guide Wire 1.0mm x 225mm	3

Implants

Code No.	Description	Qty.
T316	Herbert Screw 3.5mm Cannulated Length : 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm	4
T316A	Herbert Screw 2.4mm Cannulated Length : 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm	4

815.001 Quick Coupling Stright Handle For Headless Screw



815.002 Protection Drill Sleeve 6.5mm x 5.2mm



815.003 Guide Wire Sleeve 5.2mm x 1.5mm



815.004 Open Measuring Device For Headless Screw



815.005 Screw Holding Forcep



815.006 Cannulated Drill Bit Ø1.8mm x 90mm for Ø2.5mm Screw With Quick Coupling Handle

815.013 Cannulated Drill Bit Ø2.2mm x 115mm for Ø3.0mm Screw with Quick Coupling Handle



Hollow Conical Drill Bit With Quick Coupling Handle



Code No.	Drill Bit Diameter	Drill Bit Length	For Screw Diameter
815.007	1.8mm	90mm	2.5mm
815.014	2.2mm	115mm	3.0mm
815.019	2.5mm	115mm	3.0mm
815.023	2.8mm	115mm	4.0mm
815.030	3.2mm	115mm	5.0mm

815.008 Canulated Screw Driver 1.5mm Tip for Ø2.5mm Screw with Quick Coupling Handle

815.020 Canulated Screw Driver 2.0mm Tip for Ø3.0/3.5mm Screw with Quick Coupling Handle

815.024 Canulated Screw Driver 2.5mm Tip for Ø4.0mm Screw with Quick Coupling Handle

815.031 Canulated Screw Driver 3.0mm Tip for Ø5.0mm Screw with Quick Coupling Handle



815.009 Screw Driver 1.5mm Tip for Ø2.5mm Screw with Quick Coupling Handle

815.021 Screw Driver 2.0mm Tip for Ø3.0/3.5mm Screw with Quick Coupling Handle

815.025 Screw Driver 2.5mm Tip for Ø4.0mm Screw with Quick Coupling Handle

815.032 Screw Driver 3.0mm Tip for Ø5.0mm Screw with Quick Coupling Handle



815.010 Drill Guide 0.8mm for Ø2.5mm Screw

815.015 Drill Guide 1.2mm for Ø3.0/3.5mm Screw

815.026 Drill Guide 1.5mm for Ø4.0/5.0mm Screw



815.011 Guide Wire 0.8mm for Ø2.5mm Screw

815.016 Guide Wire 1.2mm for Ø3.0/3.5mm Screw

815.027 Guide Wire 1.5mm for Ø4.0/5.0mm Screw



815.012 Guide wire 0.8mm With Eye for Ø2.5mm Screw

815.017 Guide wire 1.2mm With Eye for Ø3.0/3.5mm Screw

815.028 Guide wire 1.5mm With Eye for Ø4.0/5.0mm Screw




815.018 Cannulated Drill Bit Ø2.5mm x 150mm for Ø3.5mm Screw

815.022 Cannulated Drill Bit Ø2.8mm x 150mm for Ø4.0mm Screw

815.029 Cannulated Drill Bit Ø3.2mm x 150mm for Ø5.0mm Screw



SAGC.26 Headless Screw Implant And Instrument Set



Implants

Code No.	Description	Qty.
T2001	Headless Screw 2.5mm Cannulated Length : 08mm, 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm	1
T2002	Headless Screw 3.0mm Cannulated Length : 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 35mm, 40mm, 45mm	1
T2003	Headless Screw 3.5mm Cannulated Length : 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 35mm, 40mm, 45mm, 50mm	1
T2004	Headless Screw 4.0mm Cannulated Length : 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 35mm, 40mm, 45mm, 50mm	1
T2005	Headless Screw 5.0mm Cannulated Length : 30mm, 35mm, 40mm, 45mm, 50mm, 55mm, 60mm, 65mm, 70mm, 75mm, 80mm	1

Code No.	Description	Qty.
815.001	Quick Coupling Stright Handle For Headless Screw	1
815.002	Protection Drill Sleeve 6.5mm x 5.2mm	1
815.003	Guide Wire Sleeve 5.2mm x 1.5mm	1
815.004	Open Measuring Device For Headless Screw	1
815.005	Screw Holding Forcep	1
815.006	Cannulated Drill Bit Ø1.8mm x 90mm for Ø2.5mm Screw With Quick Coupling Handle	2
815.007	Hollow Conical Drill Bit Ø1.8mm x90mm for Ø2.5mm Screw With Quick Coupling Handle	1
815.008	Canulated Screw Driver 1.5mm Tip for Ø2.5mm Screw with Quick Coupling Handle	1
815.009	Screw Driver 1.5mm Tip for Ø2.5mm Screw with Quick Coupling Handle	1
815.010	Drill Guide 0.8mm for Ø2.5mm Screw	1
815.011	Guide Wire 0.8mm for Ø2.5mm Screw	4
815.012	Guide wire 0.8mm With Eye for Ø2.5mm Screw	2
815.013	Cannulated Drill Bit Ø2.2mm x 115mm for Ø3.0mm Screw with Quick Coupling Handle	2
815.014	Hollow Conical Drill Bit Ø2.2mm x 115mm For Ø3.0mm Screw with Quick Coupling Handle	1
815.015	Drill Guide 1.2mm for Ø3.0/3.5mm Screw	2
815.016	Guide Wire 1.2mm for Ø3.0/3.5mm Screw	4
815.017	Guide wire 1.2mm With Eye for Ø3.0/3.5mm Screw	2
815.018	Cannulated Drill Bit Ø2.5mm x 150mm for Ø3.5mm Screw	2

Code No.	Description	Qty.
815.019	Hollow Conical Drill Bit Ø2.5mm x 115mm For Ø3.0mm Screw with Quick Coupling Handle	1
815.020	Canulated Screw Driver 2.0mm Tip for Ø3.0/3.5mm Screw with Quick Coupling Handle	2
815.021	Screw Driver 2.0mm Tip for Ø3.0/3.5mm Screw with Quick Coupling Handle	2
815.022	Cannulated Drill Bit Ø2.8mm x 150mm for Ø4.0mm Screw	2
815.023	Hollow Conical Drill Bit Ø2.8mm x 115mm For Ø4.0mm Screw with Quick Coupling Handle	1
815.024	Canulated Screw Driver 2.5mm Tip for Ø4.0mm Screw with Quick Coupling Handle	1
815.025	Screw Driver 2.5mm Tip for Ø4.0mm Screw with Quick Coupling Handle	1
815.026	Drill Guide 1.5mm for Ø4.0/5.0mm Screw	2
815.027	Guide Wire 1.5mm for Ø4.0/5.0mm Screw	4
815.028	Guide wire 1.5mm With Eye for Ø4.0/5.0mm Screw	2
815.029	Cannulated Drill Bit Ø3.2mm x 150mm for Ø5.0mm Screw	2
815.030	Hollow Conical Drill Bit Ø3.2mm x 115mm For Ø5.0mm Screw with Quick Coupling Handle	1
815.031	Canulated Screw Driver 3.0mm Tip for Ø5.0mm Screw with Quick Coupling Handle	1
815.032	Screw Driver 3.0mm Tip for Ø5.0mm Screw with Quick Coupling Handle	1



House of Ortho Tools

Hardik International Pvt. Ltd.

Plot No. G-2041/42/A, Kishan Gate, Lodhika G.I.D.C.,
Kalawad Road, Metoda, RAJKOT-360 021, (INDIA).

info@hardikinternational.com
www.hardikinternational.com

☎ : +91 97370 47470/1/2/4 | 📞 : hardikinternational

Follow Us :     