



FEMORAL NECK SYSTEM

IMPLANTS & INSTRUMENTS

T271B Femoral Neck System Plate 130°




HOLES

1hole, 2hole

MATERIAL

Titanium

T272 Femoral Neck System Bolt



LENGTH

75mm to 130mm (5mm Increment)

MATERIAL

Titanium

T273 Antirotation Screw for FNS



LENGTH

75mm to 130mm (5mm Increment)

MATERIAL

Titanium

T225 5.0mm Locking Compression Screws (Self Tapping)




LENGTH

30mm to 60mm (2mm Increment)

MATERIAL

Titanium

T466 5.0mm Locking Compression Screws (Self Tapping) Star Head



LENGTH

30mm to 60mm (2mm Increment)

MATERIAL

Titanium

814.001 Insertion Handle for FNS



814.002 Insertion Handle Device for FNS



814.003 Protection Sleeve for FNS



814.004 Multifunction Rod Insertion for FNS



814.005 Depth Gauge up to 100mm for FNS



814.006 Direct Measuring Device for 3.2mm Guide Wire for FNS



814.007 Reamer for FNS



814.008 130° Angled Guide for 3.2 mm Guide Wire for FNS



814.009 Correction Guide for 3.2 mm Guide Wire for FNS



814.010 Hammer Guide for FNS



814.011 Guide Pin Resetter



814.012 Quick Coupling Handle for FNS



735.08 Torque Limit Adapter 4.0 Nm



612.011 Quick Coupling Stright Handle



814.015 Shaft for Star Screw Driver T25 for FNS



814.016 Drill Stop for FNS



814.017 Threaded Guide Wire 3.2mm
Length 400mm for FNS



814.018 Guide Pin for FNS



814.019 Drill Bit 4.3mm x 413 mm for FNS



814.020 Shaft for Hex Screw Driver 3.5mm Tip for FNS



SAGC.20 Femoral Neck System Set



Implants

Code No.	Description	Qty.
T271B	Femoral Neck System Plate 130° Holes : 1hole, 2hole	2
T272	Femoral Neck System Bolt Length : 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm, 125mm, 130mm	1
T273	Antirotation Screw for FNS Length : 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm, 125mm, 130mm	1
T225	5.0mm Locking Compression Screws (Self Tapping) Length : 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm	2
OR		
T466	5.0mm Locking Compression Screws (Self Tapping) Star Head Length : 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm	2

Instruments

Code No.	Description	Qty.
814.001	Insertion Handle for FNS	1
814.002	Insertion Handle Device for FNS	1
814.003	Protection Sleeve for FNS	1
814.004	Multifunction Rod Insertion for FNS	1
814.005	Depth Gauge up to 100 mm for FNS	1
814.006	Direct Measuring Device for 3.2mm Guide Wire for FNS	1
814.007	Reamer for FNS	1
814.008	130° Angled Guide for 3.2 mm Guide Wire for FNS	1
814.009	Correction Guide for 3.2 mm Guide Wire for FNS	1
814.010	Hammer Guide for FNS	1
814.011	Guide Pin Resetter	1
814.012	Quick Coupling Handle for FNS	1
735.08	Torque Limit Adapter 4.0 Nm	1
612.011	Quick Coupling Stright Handle	1
814.015	Shaft for Star Screw Driver T25 for FNS	2
814.016	Drill Stop for FNS	1
814.017	Threaded Guide Wire 3.2mm Length 400mm for FNS	5
814.018	Guide Pin for FNS	2
814.019	Drill Bit 4.3mm x 413 mm for FNS	1
	Optional	
814.020	Shaft for Hex Screw Driver 3.5mm Tip for FNS	2



House of Ortho Tools

Hardik International Pvt. Ltd.

Plot No. G-2041/42/A, Kishan Gate, Lodhika G.I.D.C.,
Kalawad Road, Metoda, RAJKOT-360 021, (INDIA).

info@hardikinternational.com
www.hardikinternational.com

☎ : +91 97370 47470/1/2/4 | 📞 : hardikinternational

Follow Us :     