

Hardik International Private Limited

PRODUCT COMBINATION CHART

Doc. No.: D MK 10  
 Origin Date : 01.08.2020  
 Rev. No.: 00  
 Rev. Date : -----

System Name	Product Name	Use with main product as combination (Product or Accessory)					Instrument Set code	Instrument Set Name	Instrument Set Details	Mention if used on joints (Implants used on joints will affect the adjacent body part)
	(Plate, Nail etc.)	Screws(If Applicable)With Quantity	Pins (If applicable) With Quantity	Wire (If Applicable) With Quantity	Bolt (If Applicable)With Quantity	Applicable) With Quantity				
2.7mm Plating system(non-Locking)	2.7mm DCP Plate	Cortex Screw								
	2.7mm Quarter Tubular Plate	Self Tapping Cortex Screw								
	2.7mm L-Plate 2Hole (Left & Right)	Cortex Screw Self Tapping Star Head								
2.7mm Plating system(Locking)	Locking Volar Distal Radius Plate (Left & Right)	Cortex Screw								
	Locking XVA Volar Distal Radius Plate (Left & Right) Head 6 hole Narrow	Self Tapping Cortex Screw								
	Locking XVA Volar Distal Radius Plate (Left & Right) Head 6 hole	Cortex Screw Self Tapping Star Head								
	Locking XVA Volar Distal Radius Plate (Left & Right) Head 7 hole	Locking Compression Screw								
	Locking Proximal Radius Neck Plate	Locking Compression Screw Star Head								
	Locking Proximal Radius Rim Plate (Left & Right)	Locking Compression Screw Star Head (Variable Angle)								
	Locking Distal Radius Volar Rim Plate (Left & Right) Head 6 Hole									
	Locking Distal Radius Volar Rim Plate (Left & Right) Head 7 Hole		N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	Locking Distal Radius Plate (Left & Right)									
	Locking Distal Radius Buttress Plate 2.7mm (Head 5hole) (Left & Right)									
	Locking Distal Radius Plate Extra Articular (Head 5hole) (Left & Right)									
	Locking Distal Radius Plate 2.7mm Extra Articular (Head 4hole) (Left & Right)									
	Locking Distal Radius Plate Straight									
Locking Distal Radius L Plate (Head 2hole) (Left & Right)										

Locking Distal Radius L Plate (Head 3hole) (Left & Right)									
Locking Distal Radius L Plate Oblique Angled(Head 3hole) (Left & Right)									
Locking Distal Radius T Plate (Head 3hole)									
Locking Distal Ulna Hook Plate									
Small Locking Compression Plate	Locking Cancellous Screw (Short Thread)								
Locking Compression T Plate, Right Angled	Locking Cancellous Screw (Short Thread) Star Head								
Locking T Plate Oblique Angled (Left & Right)	Cortex Screw								
Locking Reconstruction Plate	Self Tapping Cortex Screw								
Locking Proximal Humeral Plate	Cortex Screw 14 TPI								
Locking Medial Distal Tibia Plate (Left & Right)	Cortex Screw Self Tapping Star Head								
Locking Cloverleaf Plate	Locking Compression Screw								
Locking Long Head Proximal Humeral Plate	Locking Compression Screw Star Head								
Locking Long Head Proximal Humeral Plate (Broad)	Locking Compression Screw Star Head (Variable Angle)								
Locking Distal Humeral Plate with Support 3.5mm (Left & Right)									
Locking Distal Humeral Plate (Left & Right)									
Locking Medial Distal Humeral Plate (Left & Right)									
Locking Superior Anterior Clavicle Plate (Left & Right)									
Locking Superior Anterior Clavicle Plate with Lateral Extension (Left & Right)									
Locking Clavicle Hook Plate Hook Depth (12mm to 18mm) (Left & Right)									



3.5mm Plating system(non-Locking)	Reconstruction Straight Plate (DC Hole)	Cortex Screw	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	3.5mm Reconstruction Plate (Round Hole)	Self Tapping Cortex Screw								
	Pelvic Reconstruction Plate	Cortex Screw 14 TPI								
	Straight Reconstruction Plate(Thin)	Cortex Screw Self Tapping Star Head								
	One Third Tubular Plate	Cancellous Screw Short Thread								
	Small DCP Plate	Cancellous Screw Short Thread Star Head								
	Small Limited Dynamic Compression Plate	Cancellous Screw (Full Thread)								
	Small T-Ellis Plate	Cancellous Screw (Full Thread) Star Head								
	Small T Plate	Malleolar Screw								
	Small T Plate Oblique Angle									
	Small T Plate Right Angled (Head 4 hole)									
	Small T Plate Right Angled									
	Small Hook Plate									
	Cloverleaf Plate									
	Calcaneal Plate									
Semi One-Third Compression Plate										
	Straight Reconstruction Plate(Thin)	Cortex Screw								
	T Plate	Self Tapping Cortex Screw								
	T-Plate	Cortex Screw Self Tapping Star Head								

4.5/5.0mm Plating System(Non Locking)	Small Hook Plate	Cancellous Screw (16 mm Thread)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Semi Tubular Plate	Cancellous Screw (32mm Thread)								
	Narrow DCP Plate	Cancellous Screw (Full Thread)								
	Narrow Limited Contact Dynamic Compression Plate	Cancellous Screw (Full Thread) Star Head								
	Broad Limited Contact Dynamic Compression Plate	Malleolar Screw								
	Broad DCP Plate									
	T-Buttress Plate									
	L-Buttress Plate (Left & Right)									
	L- Plate (Left & Right)									
	Spoon Plate									
	Lateral Tibia Head Buttress Plate (Left & Right)									
	Condylar Buttress Plate (Left & Right)									
Cobra Head Plate										
	Straight Reconstruction Plate(Thin)	Locking Cancellous Screw (16mm Thread)								
	T Plate	Locking Cancellous Screw (16mm Thread) Star Head								
	T-Plate	Locking Cancellous Screw (32mm Thread)								
	Small Hook Plate	Locking Cancellous Screw (32mm Thread) Star Head								
	Semi Tubular Plate	Locking Cancellous Screw (Full thread)								
	Narrow DCP Plate	Locking Cancellous Screw (Full thread) Star Head								

4.5/5.0mm Plating System(Locking)

Narrow Limited Contact Dynamic Compression Plate	Locking Compression Screw
Broad Limited Contact Dynamic Compression Plate	Locking Compression Screw Star Head
Broad DCP Plate	Locking Compression Screw Star Head (Variable Angle)
T-Buttress Plate	Cortex Screw
L-Buttress Plate (Left & Right)	Self Tapping Cortex Screw
L- Plate (Left & Right)	Cortex Screw Self Tapping Star Head
Spoon Plate	
Lateral Tibia Head Buttress Plate (Left & Right)	
Condylar Buttress Plate (Left & Right)	
Cobra Head Plate	
Locking Reconstruction Plate	
Narrow Locking Compression Plate	
Broad Locking Compression Plate	
Locking Reconstruction Plate	
Locking L-Buttress Plate(Left & Right)	
Locking T-Buttress Plate	
Locking T- Plate	
Locking Lateral Tibia Plate (Left & Right)	
Locking Proximal Tibia Plate (Left & Right)	

N/A

N/A

N/A

N/A

N/A

N/A

N/A

N/A

	Locking Distal Femoral Plate (Left & Right)										
	Locking Condylar Buttress Plate (Left & Right)										
	Locking Proximal Femoral Plate (Left & Right)										
	Locking Proximal Femoral Plate (Left & Right) Head 4Hole										
	Locking Tomo Fix Plate										
	Locking Distal Medial Femoral Plate (Left & Right)										
3.5mm/4.5mm/5.0mm Plating System(Locking)	Locking Metaphyseal Plate										
DHS PLATES NON LOCKING	DHS Plate (DC Hole) 38mm Barrel 125° to 150°	Compression For DHS/DCS Lag Screw	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	DHS Plate (DC Hole) 32mm Barrel 125° to 150°	DHS/DCS Lag Screw									
	DHS Plate (DC Hole) 25mm Barrel 125° to 140°	Pediatric DHS/DCS Lag Screw									
	DHS Plate (Round Hole) 38mm Barrel 125° to 140°	Cortex Screw									
	DHS Plate (Round Hole) 32mm Barrel 125° to 140°	Self Tapping Cortex Screw									
	DHS Plate (Round Hole) 25mm Barrel 125° to 140°	Cortex Screw Self Tapping Star Head									
	95° DCS Plate (DC Hole)										
	95° DCS Plate 19mm (DC Hole)										
	95° DCS Plate (Round Hole)										
	Pediatrics DHS Plate (DC Hole) 20mm Barrel 130° to 135°										
	Trochanter Stabilization Plate for DHS										
	Locking DHS Plate 38mm Barrel 125° to 140°	Compression For DHS/DCS Lag Screw									

DHS PLATES LOCKING	Locking DHS Plate 25mm Barrel 125° to 140°	DHS/DCS Lag Screw	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Locking 95° DCS Plate	Cortex Screw								
		Self Tapping Cortex Screw								
		Cortex Screw Self Tapping Star Head								
		Locking Compression Screw								
FOOT AND ANKLE PLATES	Eight Plate	Cortex Screw	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Eight Plate Z Angle	Self Tapping Cortex Screw								
	First MTP Fusion Plate (Left & Right)	Locking Compression Screw								
	First TMT Fusion Plate	Locking Compression Screw Star Head								
	Second TMT Fusion Plate	Cortex Screw Self Tapping Star Head								
	Cloverleaf Fusion Plate									
	T Fusion plate									
	L Fusion Plate (Left & Right)									
	Straight Fusion Plate									
	X – Plate									
	Foot Calcaneal Plate (Left & Right)									
FNS PLATES	Femoral Neck System Plate 130° to 135°	Antirotation Screw for FNS	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
		Femoral Neck System Bolt								
		Locking Compression Screw								



		Locking Compression Screw Star Head								
MAXILOFACIAL PLATES NON LOCKING	1.5mm/ 2.0mm/ 2.4mm Mandible Plate	Cortex Screw (Self Tapping) Cross Head	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	1.5mm/ 2.0mm/ 2.4mm Mandible Plate with Center Space	Cortex Screw Self Drilling Cross Head								
	1.5mm/ 2.0mm/ 2.4mm Orbital Plate with Centre Space	Emergency Screw Cross Head								
	1.5mm/ 2.0mm/ 2.4mm Orbital Plate	IMF Screw								
	1.5mm/ 2.0mm/ 2.4mm L Plate Angle 2hole (Left & Right)	Cortex Screw								
	1.5mm/ 2.0mm/ 2.4mm L Plate Angle 3hole (Left & Right)	Self Tapping Cortex Screw								
	1.5mm/ 2.0mm/ 2.4mm T Plate Head 2Hole									
	1.5mm/ 2.0mm/ 2.4mm X Plate									
	1.5mm/ 2.0mm/ 2.4mm Y Plate									
	1.5mm/ 2.0mm/ 2.4mm Curved Plate (Left & Right)									
	1.5mm/ 2.0mm/ 2.4mm L Plate Angle 2 (Without Bar) (Left & Right)									
	1.5mm/ 2.0mm/ 2.4mm Angular Plate (Left & Right)									
	1.5mm/ 2.0mm/ 2.4mm "L" Shape 120° Oblique Angled Plate (Left & Right)									
	1.5mm/ 2.0mm/ 2.4mm Nose "Y" Shape Plate									
	1.5mm/ 2.0mm/ 2.4mm Nose "H" Shape Plate									
	1.5mm/ 2.0mm/ 2.4mm Chin Plate									
	1.5mm/ 2.0mm/ 2.4mm I Plate									
	1.5mm/ 2.0mm/ 2.4mm 3D Plate									

1.5mm Orbital Floor Plate								
2.0mm / 2.4mm Double End Y Plate								
2.0mm / 2.4mm Z Plate (Left & Right)								
2.4mm/ 2.7mm Angled Reconstruction Plate								
2.4mm Reconstruction Plate Full Mandible Shape								
2.0mm Delta Plate								
Titanium Mesh 1.5/2.0/2.4mm								
Locking Mesh	Locking Cortex Screw (Self Tapping) Cross Head							
Locking 1.5mm/2.0mm/2.4mm Mandible Plate with Center Space	Locking Cortex Screw Self Drilling Cross Head							
Locking 1.5mm/2.0mm/2.4mm Mandible Plate	Locking Cortex Screw							
Locking 1.5mm/2.0mm/2.4mm Orbital Plate with Centre Space								
Locking 1.5mm/2.0mm Orbital Plate								
Locking 1.5mm/2.0mm/2.4mm L Plate Angle 2 (Left & Right)								
Locking 1.5mm/2.0mm L Plate Angle 3 (Left & Right)								
Locking 1.5mm/2.0mm/2.4mm T Plate Head 2								
Locking 1.5mm/2.0mm/2.4mm X Plate								
Locking 1.5mm/2.0mm/2.4mm Y Plate								
Locking 1.5mm/2.0mm/2.4mm 3D Plate								
Locking 1.5mm L Plate Angled 2Hole (Left & Right) (Without Bar)								

MAXILOFACIAL PLATES LOCKING	Locking 1.5mm/2.0mm/2.4mm L Plate Angled 2Hole (Left & Right)		N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Locking 1.5mm/2.0mm/2.4mm Angular Plate 150° (Left & Right)									
	Locking 1.5mm/2.0mm L Plate Angled 3Hole (Left & Right) (Without Bar)									
	Locking 1.5mm/2.0mm/2.4mm L Plate Angled 3Hole (Left & right)									
	Locking 1.5mm/2.0mm/2.4mm "L" Shape 120° oblique Angle Plate (Left & Right)									
	Locking 1.5mm/2.0mm/2.4mm Nose "Y" Shape Plate									
	Locking 1.5mm/2.0mm/2.4mm Nose "H" Shape Plate									
	Locking 1.5mm/2.0mm/2.4mm Chin Plate									
	Locking 1.5mm/2.0mm/2.4mm I Plate									
	Locking 2.0mm/2.4mm Double End Y Plate									
	Locking 2.0mm/2.4mm Z Plate (Left & Right)									
	Locking 2.0mm/2.4mm Curved Plate (Left & Right)									
	Locking 2.4mm TMJ Plate (Left & right)									
	Locking 2.4mm Angled Reconstruction Plate (Left & Right)									
Universal Interlocking Nail Tibia	End Cap for Universal Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA	
Universal Interlocking Nail Femoral	End Cap for Universal Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA	
XL Tibia Interlocking Nail	End Cap for XI Nail Tibia & Femoral	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA	
XL Interlocking Femoral Nail	End Cap for XI Nail Tibia & Femoral	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA	
XL Tibia Interlocking Nail (Distal Tip)	End Cap for XI Nail Tibia & Femoral	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA	

Nail	Supra Condylar Nail	End Cap for Supra Condylar Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Humeral Nail	End Cap for Humeral Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Multilock Humeral Nail (Left & Right)	End Cap for Multilock Humeral nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Reconstruction Femoral Nail (Left & Right)	End Cap for Reconstruction Femoral Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Proximal Femoral Nail 130° & 135°	6.2mm Hip Locking Bolt ,8.0mm Proximal Femoral Screw, End Cap for Proximal Femoral Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Long Proximal Femoral Nail (Left & Right) 130° & 135°	6.2mm Hip Locking Bolt ,8.0mm Proximal Femoral Screw, End Cap for Proximal Femoral Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Trochanteric Femoral Nail 130° & 135°	6.2mm Hip Locking Bolt ,8.0mm Proximal Femoral Screw, End Cap for Proximal Femoral Nail	NA	NA	Interlocking Bolt	NA	NA	NA	NA	NA
	Elastic Nail	End Cap For Elastic Nail	NA	NA	N/A	NA	N/A	N/A	N/A	N/A
	Enders Nail	N/A	NA	NA	N/A	N/A	N/A	N/A	N/A	N/A
	K Nail	N/A	NA	NA	N/A	N/A	AGC.102	N/A	N/A	N/A
	Rush Nail	N/A	NA	NA	N/A	N/A	N/A	N/A	N/A	N/A
	Square Nail Radius	N/A	NA	NA	N/A	N/A	N/A	N/A	N/A	N/A
	Square Nail Ulna	N/A	NA	NA	N/A	N/A	N/A	N/A	N/A	N/A
	PFN-2 Nail (Left & Right) 130° & 135°	PFN-2 Blade Screw,End Cap For PFN-2 Nail,End Cap For PFN-2 Nail Locking Blade	NA	NA	Interlocking Bolt	NA	N/A	N/A	N/A	N/A
	Expert Tibia Interlocking Nail	End Cap for Expert Tibia Nail	NA	NA	Interlocking Bolt	NA	N/A	N/A	N/A	N/A
G-PFN Nail 130° & 135°	G-PFN Screw,End Cap Outer For G-PFN Nail,End Cap Inner For G-PFN Nail,End Cap For G-PFN Nail Screw	NA	NA	Interlocking Bolt	NA	N/A	N/A	N/A	N/A	
Long G-PFN Nail (Left & Right) 130° & 135°	G-PFN Screw,End Cap Outer For G-PFN Nail,End Cap Inner For G-PFN Nail,End Cap For G-PFN Nail Screw	N/A	N/A	Interlocking Bolt	N/A	N/A	N/A	N/A	N/A	
Cannulated Screw	Cannulated Screw (Short Thread)									
	Cannulated Screw (1/3 Thread)									

	Cannulated Screw (1/2 Thread)									
	Cannulated Screw (16mm Thread)									
	Cannulated Screw (32mm Thread)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Cannulated Screw (Full Thread)									
	Locking Cannulated Screw (16mm Thread)									
	Locking Cannulated Screw (32mm Thread)									
	Locking Cannulated Screw (Full Thread)									
Herbert Screw	Herbert Screw Cannulated	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Schanz Screw(16mm to 60mm Thread Length)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Schanz Screw Tapper Thread	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Cannulated Interference Screw (ACL Screw)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	Suture Anchor Screw	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A